

Hearing Loops Are Already a Proven Assistive Listening System



Today, people with hearing loss need to hear clearly in their daily activities. Here's what a person can already do with hearing loops:

- Hear better worldwide. Hearing loops are already in multiple countries.
- Connect hearing devices like hearing aids directly to the hearing loop. No intermediary device, such as a smartphone, is necessary.
- Benefit no matter their physical ability or disability, such as low vision or limited dexterity.
- Tune in easily, no matter their age.
- Be guaranteed a connection. Hearing loops can connect to an unlimited number of people with telecoil programs in their hearing devices.
- Use almost anywhere: service counters, small or large meeting rooms, indoor or outdoor environments, entertainment venues, government offices, transportation, and more.
- Switch easily between different hearing loop locations. Hearing loops may be used in adjoining areas.
- Make the most of new hearing loop technology, designs, and equipment improvements.
- Participate via a very inexpensive component (telecoils) in their hearing instruments.
- Maintain privacy and confidentiality, e.g., a person's hearing aids are not identified or logged.
- Gain the maximum benefits from proven technology. Unbiased, third-party research has been conducted, and loops work exceptionally well for people with hearing loss.
- Know that hearing loops comply with the Americans with Disabilities Act (ADA) Standards.
- Trust that hearing loops can be designed to an international standard (IEC 60118-4) for quality sound and no latency.